



Clear
Surgical



Oplight™

The Oplight

Clear Surgical have developed an innovative enabling technology, the **Oplight**. The **Oplight** is a patented, LED, sterile, lighting cartridge which fits on to existing flat blade medical retractors to light up the surgical cavity area of an operation. This solves the existing problem of shadows, light reflection and cables currently suffered by surgeons.

This information pack will provide you with all technical information you require to assist in the evaluation of the **Oplight** technology.



Powerful “Inside Surgical Cavity” Illumination

The Latest Advancement In LED “ Inside Surgical Cavity” Illumination.
Compact, Single-Use, Cable-Free, Flat Bladed Retractor Compatible!

It's as simple as “stick and click”

Advantages

- Powerful, Long Lasting Illumination
- Extremely easy to use just “Stick and Click”
- Compatible with all flat bladed retractors
- A must for ALL “Deep Surgical Cavity” operations
- On/Off switch for additional surgeon control
- Patented technology directs light into the surgical cavity and away from surgeon's eyes
- No more shadows associated with ceiling lights and head mounted lamps
- Compact and efficient (55mmx25mm actual size)
- No capital equipment expense

Frequently Asked Questions

- What class of medical device is the Oplight?

It is a Class 1s (sterile). See SGS (UK) Ltd certificate GB19/963717 Directive 93/42/EEC on medical devices, Annex V.

- What is the shelf life of the Oplight?

The Oplight has a shelf life of 3 years, as determined by the sterility process and packaging. We will always supply items for sale with at least 24 months shelf life. Samples may have a shorter shelf life.

- What size is the Oplight?

The Oplight measures 55.96mm by 27.3mm by 4.7 mm and weights 26.6g

- How will the Oplight be packed?

The Oplight is packed in single units. The Oplight will be presented in a sterile double pouched pack

- Can I switch the Oplight off and on during a procedure?

Yes. The Oplight has a single off/on switch that is accessible during use.

- Can the Oplight be taken from one retractor and placed on another during a procedure?

No. We do not advise removing the Oplight from one retractor as this may reduce the adhesion properties.

- Does the Oplight contain any latex?

The Oplight has been manufactured to be free of latex

Oplight compatible retractors

Below is a list of surgical retractors which will accommodate the Oplight during surgical interventions. Whilst this is a comprehensive list there may be additional “flat blade” retractors that will be compatible with the Oplight.

Langenbeck Retractor	Kocher Retractor
Phyris Retractor	Tudor Edwards Scapula Retractor
St Marks Retractor	Cloward Style Retractor
Norfolk and Norwich Retractor	Bookler Right Angle Retractor
Obwegeser Retractor	Jackson Burrows Retractor
Simon Vaginal Retractor OB/GYN	Hibbs Retractor
Pryor-Pean Retractor OB/GYN	Taylor Spinal Retractor
Doyen Vaginal Retractor OB/GYN	Meyerding Retractor
Heaney Hysterectomy Retractor OB/GYN	Lamina Spreader
Guttman OB/GYN	Fakuda Shoulder Retractor
Lateral Vaginal Wall Retractor OB/GYN	Tuffier Retractor
Reinhoff Retractor	Lemmon Extra Deep Deaver Retractor
Heaney-Simon Retractor	Yu Holtgreve Malleable Blade
Piling Multipurpose blade	Blafour Center Blades
Davidson Scapula Retractor	Wexler Expandable Blade
Harrington Extra Deep	Beckman-Eaton Retractor
Deep Deaver Retractor	Davis Brain Spatula
DeBakey- Cooley Retractor	Sauerbruch Retractor
Taylor Retractor	Wishbone Splachnic Retractor
Horizontal Retractor	

Oplight Surgeries

The Oplight has been used in number of surgeries which have been included below. This list is not exhaustive.

Vascular

Adominal aortic endarterectomy
Fem-pop bypass in a deep popliteal fossa
Nerogenic thoracic outlet syndrome
Vascular thoracic outlet syndrome
Popliteal aneurism

ENT

Parotidectomy
Thyroidglossal Cysts
Branchial Cysts
Tonsillectomy

Colorectal

Anal fissures
Anal fistula

Others

Breast & axillary access
Axillary node dissection
Mastectomy; depending on approach used
Plastic Surgery
Cardiac TAVI surgery
Deep posterior capsule release
Anterior Lumbar Interbody Fusion (ALIF)
Lumbar posterior fusion (PLIF)
Open Cholecystectomy
Open Appendectomy
Hernia
Carpal Tunnel
Anterior Hip Arthroplasty

Oplight Specification Rev 1.1

Disposable single-use, sterile cartridge incorporating LED technology to provide light source to directly illuminate the interior cavity of an incision during surgery.

Power Source:

Li-MnO₂ Battery - CP162235/160mAh/3.0V

Material:

Resin - ALTUGLAS SG7

Adhesive – Avery Dennison MED 6361U

Light Source

NSSW206 White LED

Storage


Storage Temperature Range: 15°C to 27°C

Storage Humidity Range: 30% to 60% relative humidity (RH)

Atmospheric pressure: 97kPa to 104 kPa

Sizes and Codes

SLC5525DT - Sterile



Product Weight

Light Cartridge 26.6g

Performance Characteristics

The average illumination over 45mins - 1132 Lux.

Caution

- Ensure there is no damage to sterile barrier packaging.
- Do not use if packaging is damaged.
- Only to be used by appropriately trained individuals.
- Do not misuse.
- Do not use after use by date.
- Single use only.
- Do not re-sterilise.
- Dispose of Cartridge as a contaminated medical device containing a lithium battery.
- No modification of this equipment is allowed.

Guidance and Manufacturer's declaration – electromagnetic emissions

Oplight is intended for use in the electromagnetic environment specified below. The user of Oplight should assure that it is used in such an environment

Emissions Test	Compliance	Electromagnetic Environment - guidance
RF emissions CISPR 11	Group 1	Oplight uses RF energy only for its internal function. Therefore its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic Emissions IEC6100-3-2	Not Applicable	
Voltage fluctuations/flicker emissions IEC6100-3-2	Not Applicable	

Guidance and Manufacturer's declaration – electromagnetic immunity

Oplight is intended for use in the electromagnetic environment specified below. The user of Oplight should assure that it is used in such an environment

Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	+/- 6kV contact +/- 8kV air	+/- 6kV contact +/- 8kV air	Floors should be wood, concrete or ceramic. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst	+/- 2kV for power supply lines	Not Applicable	

IEC 61000-4-4	+/- 1kV for input/output lines		
Surge IEC 61000-4-5	+/- 1kV line(s) to line(s) +/- 2kV line(s) to earth	Not Applicable	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle <5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 0.5 s	Not Applicable	
Power Frequency(50/60Hz) Magnetic Field IEC 61000-4-8	3A/m	3A/m 50 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial hospital or environment
NOTE U_T is the a.c. mains voltage prior to application of the test level			

The device was tested under Classification 4 Group 1 Class B of EN 60601-1-2:2007 and was tested to the following standards:

EN 60601-1-2:2007
EN 60601-2-41: 2009

Emissions

EN 55011 Class B Radiated Disturbance

Immunity

EN61000-4-3 Radiated RF electromagnetic fields
EN61000-4-2 Electrostatic Discharge (ESD)
EN61000-4-8 Power frequency magnetic fields

Iver Macdonald

Technical & Regulatory Manager

Date: 18th October 2019

Oplight EC Declaration of Conformity Issue 1.7

This is to confirm that the Oplight is classified as a Class Is medical device in accordance with Annex IX Rule 7 and Rule 12 of the Medical Device Directive 93/42/EEC of June 1993 and amending Directive 2007/47/EEC. The Oplight complies with the essential requirements as detailed in Annex I of the Directive and is CE marked in accordance with Annex II of the Directive.

EC Certificate Full Quality Assurance System: Certificate GB19/963717

The product will be marketed in the European Economic Area (EEA) as:

Oplight

The manufacturer responsible for placing the Oplight on the Market is Clear Surgical Ltd, Glenbervie Business Centre, Ramoyle House, Glenbervie Business Park, Larbert, FK5 4RB, United Kingdom.

The route to conformity has been assessed by SGS United Kingdom Ltd, Unit 202b, Worle Parkway, Weston-super-Mare, BS22 6WA, United Kingdom.

Notified Body Number 0120

Signed:

Dated: 11/10/2019



Michelle Ward
Chairman & CEO

Clear Surgical Ltd. Glenbervie Business Centre, Ramoyle House, Glenbervie Business Park
Larbert, FK5 4RB, United Kingdom

BS EN 556-1:2001

Sterilization of medical devices. Requirements for medical devices to be designated "STERILE".
Requirements for terminally sterilized medical devices.

BS EN 1041:2008+A1:2013

Information supplied by the manufacturer of medical devices.

BS EN ISO 10993-1:OCTOBER 2009

Biological evaluation of medical devices –
Part 1: Evaluation and testing (ISO 10993-1:2003)

BS EN ISO 11135:2014

Sterilization of health care products. Ethylene oxide. Requirements for development, validation and routine control of a sterilization process for medical devices.

BS EN ISO 11607-1:2009+A1:2014

Packaging for terminally sterilized medical devices. Requirements for materials, sterile barrier systems and packaging.

BS EN ISO 11737-1:2006

Sterilization of medical devices. Microbiological methods. Determination of a population of microorganisms on products.

BS EN ISO 13485:2016

Medical devices – Quality management systems – Requirements for regulatory purposes.

BS EN ISO 14155:2011

Clinical investigation of medical devices for human subjects. Good clinical practice.

BS EN ISO 14937:2009

Sterilization of health care products. General requirements for characterization of a sterile agent and the development, validation and routine control of sterilization process for medical devices.

BS EN ISO 14971:2012

Medical devices – application of risk management to medical devices.

BS EN ISO 15223-1:2016

Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. General requirements.

BS EN ISO 15883-1:2009+A1:2014

Washer – disinfectors. General requirements, terms and definitions and tests.

BS EN 60601-1:2006+A12:2014

Medical electrical equipment. General requirements for basic safety and essential performance

BS EN 60601-1-2:2007

Medical electrical equipment. General requirements for basic safety and essential performance. Collateral standard. Electromagnetic compatibility. Requirements and tests.

BS EN 60601-2-41:2009

Medical electrical equipment. Particular requirements for basic safety and essential performance of surgical luminaires and luminaires for diagnosis.

BS EN 55011:2009+A1:2010

Industrial, scientific and medical equipment. Radio-frequency disturbance characteristics. Limits and methods of measurement.

BS EN 61000-4-3:2006+A2:2010

Electromagnetic compatibility (EMC). Testing and measurement techniques. Radiated, radio-frequency, electromagnetic field immunity test.

BS EN 61000-4-2:2009

Electromagnetic compatibility (EMC). Testing and measurement techniques. Electrostatic discharge immunity test.

BS EN 61000-4-8:2010

Electromagnetic compatibility (EMC). Testing and measurement techniques. Power frequency magnetic field immunity test.

BS EN 62366-1:2015

Medical devices. Application of usability engineering to medical devices.

BS EN ISO 17664:2004

Sterilization of medical devices. Information to be provided by the manufacturer for the processing of resterilizable medical devices.

Clear Surgical ISO 13485 Certificate

Certificate GB14/92322

SGS

The management system of

Clear Surgical Limited

Glenbervie Business Centre, Ramoyle House,
Glenbervie Business Park, Larbert, FK5 4RB, UK

has been assessed and certified as meeting the requirements of

ISO 13485:2016
EN ISO 13485:2016



For the following activities

Design, Development and Manufacture of sterile disposable surgical instruments and non sterile reusable surgical instruments.

This certificate is valid from 05 December 2018 until 04 December 2020
and remains valid subject to satisfactory surveillance audits.

Re certification audit due before 30 October 2020

Issue 3. Certified since 04 December 2014

Authorised by


SGS United Kingdom Ltd
Rossmore Business Park, Ellesmere Port, Cheshire CH65 3EN UK
t +44 (0)151 350-6666 f +44 (0)151 350-6600 www.sgs.com

HC SGS 13485 2016 0118

Page 1 of 1



0005



This document is issued by the Company subject to its General Conditions of Certification Services accessible at www.sgs.com/terms_and_conditions.htm. Attention is drawn to the limitations of liability, indemnification and jurisdictional issues established therein. The authenticity of this document may be verified at

Oplight Directive 93/42/EEC certificate

CE Mark

EC Certificate Production Quality Assurance System: Certificate GB19/963717

SGS

The management system of

Clear Surgical Limited
Glenbervie Business Centre, Ramoyle House,
Glenbervie Business Park, Larbert, FK5 4RB, UK

has been assessed and certified as meeting the requirements of

Directive 93/42/EEC
on medical devices, Annex V
Restricted to the aspects of manufacture concerned with securing and
maintaining sterile conditions.

For the following products

Oplight Sterile Single Use light source

This certificate is valid from 18 July 2019 until 04 December 2022
and remains valid subject to satisfactory surveillance audits.
Re certification audit due before 30 October 2020
Issue 1. Certified since 18 July 2019

Certification is based on reports numbered GB/PC 233181

Authorised by

Jonathan M. Hall

SGS United Kingdom Ltd, Notified Body 0120
202B Worle Parkway, Weston-super-Mare, BS22 6WA, UK
t +44 (0)1934 522917 f +44 (0)1934 522137 www.sgs.com

SGS CE 18 0811

Page 1 of 1





This document is a Web version of SGS certificate for electronic use exclusively. It shall only be available by clicking on SGS Certification Mark which has been posted on Your website. It shall not be printed in anyway. This document is copyright protected. No content or appearance may be reproduced without the express written permission of SGS. Any misuse, alteration, forgery or falsification is unlawful.